

VICAT APPARATUS

VICAT APPARATUS EO/55

For determining the consistency and setting times of a normal cement mix.

Model EO/55 with set of accessories in compliance with EN European specifications.

Can also be used in other industrial fields (tableware, sanitary ware, various ceramic bodies, etc.) for various non-standardised tests.

Comprises a steel base, light alloy support, stainless steel sliding rod, adjustable index graduated scale, and counter weight.

A needle spindle or a consistency probe is fixed to the rod.

Complete with a set of 6 \varnothing 1.13 mm needles, \varnothing 10mm consistency probe, mould supporting disk, 3 plastic truncated-cone moulds with inner diameters 70x80 mm, height 40 mm and glass thermometer, scale - 10+50 °C.

Total weight 5 kg.



Art. 06CI2595

Product Code: 06CI2595

Spare Parts: Set of 6 needles –06CI2595/A

Consistency Probe – 01CI2595/C

Set of 3 truncated moulds – 06CI2595/D

Glass thermometer - 01CI2595/E